

Fibre Optic Cables and Accessories – Call 0116 2342929

PROCEDURE FOR ORDERING FIBRE OPTIC CABLE COMPLETE SETS

- 1 Choose length and core size from list below
- 2 Choose type A fitting (light Source end)
- 3 Choose type B fitting (Instrument end)
- 4 Enter all selected H numbers on your order



H595

**Zeiss Operating
Microscope
OPMI**

Special cable with two short ends, complete with two of H564 Type A Adaptors
4.8mm core 1800mm long

LIGHT SOURCE END
TYPE A

CABLE

INSTRUMENT END
TYPE B

H572

ACM & British
Eschmann
Downs



H567

GU & Storz



H568

Heine



H570

Winter & Ibe



H566

Fuji
Olympus
Rimmer



H565

Dyonics
Stryker
Wolff



H569

Pentax



3.5mm core
1800mm long

H576

4.8mm core
1800mm long

H577

3.5mm core
2300mm long

H578

4.8mm core
2300mm long

H579

H573

ACM & British



H574

Wolff



H564

Zeiss OPMI



H559

Downs



H575

Aesculap
Dyonics
Olympus, Storz, GU



H560

ACM & British
With clip



H557

Wolff with clip

